

Brain Atlas

The BRAPH 2 Developers

June 12, 2023

This is the tutorial to work with the Graphical User Interface (GUI) to manage brain atlases, which is typically the first step that is required to perform a graph analysis in BRAPH 2.0. In this Tutorial, we will explain you how to upload a brain atlas and how to visualize it.

Contents

Open the GUI	2
Upload the Brain Atlas	2
Ready Brain Atlases	3
Create a New Brain Atlas	3
Plot the Brain Atlas	4
Export the Figure	5




Figure 1: **Brain regions.** Image with a brain surface and some brain regions.




Figure 2: **Brain Atlas GUI.** Graphical user interface to work with brain atlases. *Add image with full interface.*

Open the GUI

The Brain Atlas GUI is typically the first step in the BRAPH 2.0 pipelines. You can also do it by typing `braph2`, which will open the Graphical User Interface of the BRAPH 2.0 software. Here you must first select a Pipeline containing the analyses steps that you want to apply to your data. Once a Pipeline has been selected, the first window will allow you to upload the Brain Atlas, as shown in Figure ??.

To open the Graphical User Interface and Upload the Brain Atlas, you can also do it from the command line by typing the commands in Code 1.

Code 1: Code to launch the Brain Atlas GUI. This code can be used in the MatLab command line to launch the Brain Atlas GUI.

```
1 ba = BrainAtlas(); ①
2
3 gui = GUIElement('PE', ba); ②
4 gui.get('DRAW') ③
5 gui.get('SHOW') ④
```

- ① creates a new object BrainAtlas.
- ② creates a GUI to upload the brain atlas.
- ③ draws the GUI.
- ④ shows the GUI.



Figure 3: **Brain Atlas GUI.** Graphical user interface to work with brain atlases.

Upload the Brain Atlas

In this window you have a Menu that you can use to Open a File (File -> Open) and select a Brain Atlas that is already prepared or Import (Import -> File) a Brain Atlas you have created in excel or text format. If you Import a Brain Atlas you have created and an error message appears, check again your Brain Atlas file for missing or incorrect information based on the format for new Brain Atlases we describe above. Once you have uploaded a Brain Atlas, you will notice that the information regarding the Brain Atlas ID, Brain Atlas NAME and Brain Regions sections are now filled with the information of the atlas. You can edit the information included in any of these sections.




Figure 4: **Brain Atlas GUI**. Graphical user interface to work with brain atlases.

Ready Brain Atlases

Currently, we provide several brain atlases that are commonly used in the field of brain connectomics, which can be downloaded from our website (<http://braph.org/software/brain-atlases/>).

- **AAL90** Automated Anatomical Labelling atlas with 90 cortical and subcortical regions.
- **AAL116** Automated Anatomical Labelling atlas with 116 cortical and subcortical regions, including cerebellar areas.
- **BNA** Brainnetome atlas with 246 cortical and subcortical regions.
- **Craddock** Functional atlas with 200 cortical and subcortical regions.
- **Desikan** Anatomical atlas with 68 cortical and 14 subcortical gray matter regions derived from the FreeSurfer software.
- **Destrieux** Anatomical atlas with 148 cortical and 14 subcortical gray matter regions derived from the FreeSurfer software.
- **Schaefer** Functional brain atlas with 200 cortical regions that belong to 7 different resting-state fMRI networks.

Create a New Brain Atlas

To prepare a Brain Atlas in BRAPH 2.0 format, you should create a new excel file (.xls or .xlsx), as shown in Figure 5.

Brain Atlas ID	Brain Atlas Label	Brain Atlas Notes	Brain Surface Name	X	Y	Z	Hemisphere
1	AA190						
2	AA190 Label						
3	AA190 Notes						
4	BrainMesh_ICBM152.nv						
5	PreCG.L		Precentral gyrus Left	-38.65	-5.68	50.94	left
6	PreCG.R		Precentral gyrus Right	41.37	-4.21	52.09	right
7	SFGdor.L		Superior frontal gyrus (dorsolateral) Left	-18.45	34.81	42.2	left
8	SFGdor.R		Superior frontal gyrus (dorsolateral) Right	21.9	31.12	43.82	right
9	ORBsup.L		Middle frontal gyrus (orbital part) Left	-16.56	47.32	-13.31	left
10	ORBsup.R		Middle frontal gyrus (orbital part) Right	18.49	46.1	-14.02	right
11	MFG.L		Middle frontal gyrus Left	-33.43	32.73	35.46	left
12	MFG.R		Middle frontal gyrus Right	37.59	33.06	34.04	right
13	ORBmid.L		Middle frontal gyrus (orbital part) Left	-30.65	50.43	-6.62	left
14	ORBmid.R		Middle frontal gyrus (orbital part) Right	33.18	52.59	-10.75	right
15	IFGoperc.L		Inferior frontal gyrus (opercular part) Left	-48.43	12.73	19.02	left
16	IFGoperc.R		Inferior frontal gyrus (opercular part) Right	50.2	14.98	21.41	right
17	IFGtriang.L		Inferior frontal gyrus (triangular part) Left	-45.58	29.91	13.99	left
18	IFGtriang.R		Inferior frontal gyrus (triangular part) Right	50.33	30.16	14.17	right
19	ORBinf.L		Inferior frontal gyrus (orbital part) Left	-35.98	30.71	-12.11	left
20	ORBinf.R		Inferior frontal gyrus (orbital part) Right	41.22	32.23	-11.91	right
21	ROL.L		Rolandic operculum Left	-47.16	-8.48	13.95	left
22	ROL.R		Rolandic operculum Right	52.65	-4.25	14.63	right
23	SMAL		Supplementary motor area Left	-5.32	4.85	61.38	left
24	SMAR		Supplementary motor area Right	8.62	0.17	61.85	right
25	OLF.L		Olfactory cortex Left	8.06	15.05	-11.46	left
26	OLF.R		Olfactory cortex Right	10.43	15.91	-11.26	right
27	SFGmed.L		Superior frontal gyrus (medial) Left	-4.8	49.17	30.89	left
28	SFGmed.R		Superior frontal gyrus (medial) Right	9.1	50.84	30.22	right
29	ORBsupmed.L		Superior frontal gyrus (medial orbital) Left	-5.17	54.06	-7.4	left
30	ORBsupmed.R		Superior frontal gyrus (medial orbital) Right	8.16	51.67	-7.13	right
31	RECL		Rectus gyrus Left	-5.08	37.07	-18.14	left
32	RECR		Rectus gyrus Right	8.35	35.64	-18.04	right
33	INS.L		Insula Left	-35.13	6.65	5.44	left
34	INS.R		Insula Right	39.02	6.25	2.08	right
35	ACCL		Anterior cingulate and paracingulate gyri Left	-4.04	33.4	13.95	left
36	ACCR		Anterior cingulate and paracingulate gyri Right	8.46	37.01	15.84	right

Figure 5: **Brain Atlas GUI**. Graphical user interface to work with brain atlases.

Start by writing the following information in the first 4 rows:

- Brain Atlas ID (row 1, column 1). For example: Desikan FreeSurfer
- Brain Atlas LABEL (row 2, column 1). For example: Desikan FreeSurfer Labels
- Brain Atlas NOTES (row 3, column 1). For example: Desikan FreeSurfer Nodes
- Brain Surface Name (row 4, column 1). For example: Brain-MeshICBM152.nv

Then, from row 5, you should include the IDs of the regions of your atlas (1st column), the Labels of the regions of your atlas (2nd column), the X, Y and Z coordinates (3rd, 4th and 5th columns) and the brain hemisphere or any notes you would like to add (6th column).

Plot the Brain Atlas

Once you are satisfied you can proceed with the option "Plot Brain Atlas", which will open a brain surface that contains the nodes corresponding to the regions of your atlas.

This new window has a large menu that allows you to change its visualization. We suggest you try the different options to understand how they change the figure. Importantly, within this menu there is one option called Settings Brain Surface which, when selected, will open another window, as can be seen in the snapshot below:

This window allows you to change different visualization options, which are important to create a final figure with all the nodes included in your analysis, which is often included within the 1st Figure of a manuscript.

Most things in this panel are intuitive and again we suggest that you try different options until you achieve the visualization you desire, for example:

Some things that might not be intuitive is the difference between spheres and symbols (the first one is the geometrical structure of a node, whereas the second is just a dot inside the sphere that denotes the presence of a region).

If you wish to change the size of the spheres of all nodes, you need




Figure 6: **Brain Atlas GUI.** Graphical user interface to work with brain atlases.

to change the size of one node, select other nodes in the first column, right click to select Apply to Selection.

Finally the codes for different colours in the FACECOLOR column correspond to the hexadecimal form of RGB colors, which can be found online.




Figure 7: **Brain Atlas GUI.** Graphical user interface to work with brain atlases.




Figure 8: **Brain Atlas GUI.** Graphical user interface to work with brain atlases.

Export the Figure

xxx